



CFDA CE

**BEIEN**

Focus on Minimally Invasive Surgery

## BEIEN-SURGERY™

Disposable  
Circular Staplers



### Contact us

- +0086-510-8378-2500
- +0086-510-8378-2511
- Service@Beien-Surgery.com

**Wuxi Beien Surgery Device Co.,Ltd.**

*Quality First  
Deliver life-affirming  
Result*



[www.beien-surgery.com](http://www.beien-surgery.com)

# PRODUCT INFORMATION

## INDICATIONS

Disposable circular staplers is intended to have applications throughout the alimentary tract for end-to-end, end-to-side, and side-to-side anastomoses.



## FEATURES



1. Super-Thin anvil, more convenient for procedure.
2. Bigger inner cutting diameter, reduce anastomosis stricture.
3. Rational arrangement of staples to reduce bleeding after procedure.
4. Tissue gap between anvil and cartridge, avoid tissue damage.
5. High quality titanium staples, bigger tension-resistance on anastomosis

## SPECIFICATION

Code	Description	Color	Inner DIA	Out DIA	Staple Heights
CSS17	Circular Stapler - 17mm	White	10.0mm	17.5mm	4.8mm
CSS21	Circular Stapler - 21mm	Tan	13.0mm	21.5mm	4.8mm
CSS23	Circular Stapler - 23mm	Yellow	14.5mm	23.5mm	4.8mm
CSS26	Circular Stapler - 26mm	White	16.5mm	26.0mm	4.8mm
CSS29	Circular Stapler - 29mm	Blue	19.5mm	29.0mm	4.8mm
CSS32	Circular Stapler - 32mm	Purple	22.5mm	32.0mm	4.8mm
CSS34	Circular Stapler - 34mm	Green	24.5mm	34.0mm	4.8mm