




CFDA CE

BEIEN

Focus on Minimally Invasive Surgery

BEIEN-SURGERY™

Disposable
Linear Cutter Staplers



Contact us

☎ +0086-510-8378-2500

☎ +0086-510-8378-2511

✉ Service@Beien-Surgery.com

Wuxi Beien Surgery Device Co.,Ltd.

*Quality First
Deliver life-affirming
Result*




www.beien-surgery.com

PRODUCT INFORMATION

FEATURES

1. Single knife in cartridge for each firing
2. Tissue positioning systems, ensure the tissues could not slide or overflow.
3. Balance leverage control systems, ensures parallel closure of the instrument's forks.
 - Tissue positioning Point
 - Double balancing Points
 - Thickened Steel
4. More wider stapling lip, 1.8mm/1.9mm, ensure the consistency of tissue stapling.



DISPOSABLE LINEAR CUTTER STAPLERS SPECIFICATION

Code	Description	Color	Tissue Type	Staple Heights
BLC60	Linear Cutter Staplers - 60mm	Blue	Medium	3.8mm
GLC60	Linear Cutter Staplers - 60mm	Green	Thick	4.5mm
BLC80	Linear Cutter Staplers - 80mm	Blue	Medium	3.8mm
GLC80	Linear Cutter Staplers - 80mm	Green	Thick	4.5mm
BLC100	Linear Cutter Staplers - 100mm	Blue	Medium	3.8mm
GLC100	Linear Cutter Staplers - 100mm	Green	Thick	4.5mm

CARTRIDGE UNITS SPECIFICATION

Code	Description	Color	Tissue Type	Staple Heights
BCA60	Cartridge Unit - 60mm	Blue	Medium	3.8mm
GCA60	Cartridge Unit - 60mm	Green	Thick	4.5mm
BCA80	Cartridge Unit - 80mm	Blue	Medium	3.8mm
GCA80	Cartridge Unit - 80mm	Green	Thick	4.5mm
BCA100	Cartridge Unit - 100mm	Blue	Medium	3.8mm
GCA100	Cartridge Unit - 100mm	Green	Thick	4.5mm