



CFDA CE

BEIEN

Focus on Minimally Invasive Surgery

ENDO-ARTISAN™



Disposable
Endoscopic Reload Unit
Articulating II

Contact us

+0086-510-8378-2500

+0086-510-8378-2511

Service@Beien-Surgery.com

Wuxi Beien Surgery Device Co.,Ltd.

*Quality First
Deliver life-affirming
Result*



www.beien-surgery.com

PRODUCT INFORMATION

INDICATIONS

Enabling surgeons to confidently and consistently staple in a broader range of tissue thicknesses and simplifying the cartridge selection process.



FEATURES

1. Improved burst pressure strength.
2. Reduced tissue compression stress on outer staple row.
3. Consistent hemostatic performance across a broader indicated tissue range.
4. Variable staple heights and a stepped cartridge face increase the compressibility of tissue.
5. A strong fixed anvil and target staple buckets ensure precise contact with the cartridge.
6. Provides a smooth firing force through tissue due to stronger knife bar.
7. Improved tissue retention during manipulation and transection.
8. Thick tissue purple reload can be introduced through a 12mm port.
9. Extra-Thick tissue black reload can be introduced through a 15mm port.



SPECIFICATION

Model	Description	Cartridge Color	Staple Size (Inner to outer row)
TS60P	60mm Articulating Thick Reload Unit	Purple	3.0mm, 3.5mm, 4.0mm
TS60D	60mm Articulating Extra-Thick Reload Unit	Black	4.0mm, 4.5mm, 5.0mm