



CFDA CE

BEIEN

Focus on Minimally Invasive Surgery

ENDO-ARTISAN™



Disposable
Endoscopic Stapler

Contact us

+0086-510-8378-2500

+0086-510-8378-2511

Service@Beien-Surgery.com

Wuxi Beien Surgery Device Co.,Ltd.

*Quality First
Deliver life-affirming
Result*



www.beien-surgery.com

PRODUCT INFORMATION

INDICATIONS

Combining ergonomic design, precise articulation and one-hand grasping for increased versatility - enhancing the user experience



FEATURES



1. Ten points of articulation enhance precision of approach to target tissue/vessel.
2. Surgeon-controlled articulation of up to 45 degree on either side provides improved access to surgical sites.
3. Forward and backward compatibility - the stapler can be used with all existing reload units.
4. Three shaft lengths deliver greater flexibility for procedure-specific use.

SPECIFICATION

Model	Description	Shaft Length	Order Code
ELS12-S	Disposable Endoscopic stapler -Short	60mm	H60
ELS12-C	Disposable Endoscopic stapler -Conventional	160mm	H160
ELS12-E	Disposable Endoscopic stapler -Extra	260mm	H260