



CFDA CE

BEIEN

Focus on Minimally Invasive Surgery

ENDO-ARTISAN™

Disposable
Hemorrhoidal Cutter Stapler



Contact us

+0086-510-8378-2500

+0086-510-8378-2511

Service@Beien-Surgery.com

Wuxi Beien Surgery Device Co.,Ltd.

*Quality First
Deliver life-affirming
Result*




www.beien-surgery.com

PRODUCT INFORMATION

INDICATIONS

Disposable Hemorrhoidal cutter stapler is intended to have applications throughout the anal canal to perform surgical treatment of hemorrhoid disease.



FEATURES

1. Rational arrangement of staples to reduce bleeding after procedure.
2. Cutting tissue thickness adjustable according to prolapse situation.
3. Big volume of the cartridge to make sure the enough space.
4. Tissue gap between anvil and cartridge, avoid tissue damage.

SPECIFICATION

Model	Out Diameter	Inner Diameter	Cartridge Color	Stapler Size
PPH32	32	23	Purple	4.0
PPH34	34	25	Green	4.0