

THE ASSEMBLY TYPE









Key Features :







- ◉ Shaft and handle can be disassembled.
- ◉ Disposable shaft and reusable handle.
- ◉ The shaft comes with different types of handles.
- ◉ Consistent and precise performance.
- ◉ Easy and fast assembly & disassembly.

Product Name	Picture	Product Code	Product Discription
Blunt Nose Forceps		HD1010101	laparoscopic grasper, 5mm*330mm, ratchet/without ratchet, disposable shaft
Grasping Forceps (3 teeth)		HD1010102	laparoscopic grasper, 5mm*330mm, ratchet/without ratchet, disposable shaft



Endo Clinch Grasper		HD1010103	laparoscopic grasper, 5mm*330mm, ratchet/without ratchet, disposable shaft
Fenestrated Babcock forceps		HD1010104	laparoscopic grasper, 5mm*330mm, ratchet/without ratchet, mainly used for grasping gall bladder in cholecystectomy, disposable shaft
Grasping Forceps (Serrated)		HD1010105	laparoscopic grasper, 5mm*330mm, ratchet/without ratchet, mainly used for grasping uterine in gynecological surgery, disposable shaft
Dolphin Nose Grasper		HD1010201	laparoscopic dissector, 5mm*330mm, ratchet/without ratchet, mainly used for dissecting and grasping mini tissues, disposable shaft
Mixer		HD1010202	laparoscopic dissector, right angle, 5mm*330mm, ratchet/without ratchet, disposable shaft
Maryland Grasper (30°Curved)		HD1010203	laparoscopic dissector, 30 ° curved, 5mm*330mm, ratchet/without ratchet, disposable shaft



<p>Atraumatic Intestinal Forceps</p>		<p>HD1010301</p>	<p>Laparoscopic grasper, atraumatic 5mm*330mm, ratchet/without ratchet, mainly used for grasping small intestine in gastrointestinal surgery, disposable shaft</p>
<p>Atraumatic Lung Lobe Forceps</p>		<p>HD1010302</p>	<p>laparoscopic grasper, atraumatic 5mm*330mm, ratchet/without ratchet, mainly used for grasping lung lobe, disposable shaft</p>
<p>Atraumatic Grasper</p>		<p>HD1010303</p>	<p>laparoscopic grasper, atraumatic 5mm*330mm, ratchet/without ratchet mainly used for grasping general tissues, such as liver lobe disposable shaft</p>
<p>Curved scissors</p>		<p>HD1010502</p>	<p>laparoscopic scissors, curved 5mm*330mm, ratchet/without ratchet mainly used for cutting tissues and sutures, disposable shaft</p>
<p>Handle</p>		<p>HD1010601</p>	<p>reusable handle without ratchet</p>
<p>Handle</p>		<p>HD1030601</p>	<p>monopolar laparoscopic reusable handle, without ratchet</p>



Handle without ratchet		HD1040601	reusable handle without ratchet, with silicone rotor
Handle with Ratchet		HD1020601	reusable handle with ratchet