

# STERILE DECANTER

A UNIQUE AND INNOVATIVE DESIGN



Asia  
Goldman co.  
Medical Equipment

# MALLEABILITY MAKES THE DIFFERENCE



Ideal for Use In Every Angle and Position




This simple to use product is the most economical solution for sterile fluid transfer and is manufactured exclusively in Iran by Asia Goldman

Used for Dispensing Saline and Other Fluids Safely in Clinical Environments Where Sterility is Paramount



The decanter's cap can be resealed to lead the flow of liquid. Spiked end easily puncturing most bag outlets.

50 per case .



Reference number	A	B	C	D	E
AGD01	180mm	145mm	32mm	7mm	5m
AGD02	130mm	95mm	32mm	7mm	5m